



aidapt
AUSTRALIA

Bedroom

Houghton Community Bed

The Houghton 5 function community bed is an electronically actuated acute care bed which is appropriate for adult patient care allowing, if required, the patient to adjust their own position in the bed. The bed can be fitted with side rails, IV pole and lifting pole. Powered by mains supply with secondary 24V DC battery back up system. Safe working load is 235kg and maximum patient weight is 200kg. One year warranty on parts. Three years warranty on mechanism. Back rest tilting 0-73°


- * 4 section metal mattress support
- * 11 function handset with lock function
- * Emergency stop button
- * CPR release
- * Trendelenburg & reverse trendelenburg
- * Back rest tilting 0-73°
- * Knee rest tilting 0-37°
- * 4 braked castors
- * Battery back up
- * Easy assembly
- * Brackets for transportation
- * Weight Capacity 200kg

Options: **VG891C** Metal Side Rails –pair
VG892J Lifting Pole (triangle with strap)
VG891V IV Pole
VG890B Brown or **VG890W** Light Oak end boards
 Mattress sold separately



Product Code: VG890*

Specifications



aidapt

1300 313 357
www.aidapt.com.au
shop@aidapt.com.au